

All-in-one solution for single gene disorders with diverse mutational spectrum

Thalassemia

Background

Thalassemia is a group of inherited hemolytic disorders characterized by reduced or absent production of globin chains, involving a spectrum of globin gene defects (SNVs, InDels, recombination, CNVs, etc.). It significantly impacts regions in Southeast Asia, the Mediterranean, and the Middle East.

Thalassemia carries a **high risk of missed cases** with routine blood tests

Sequential hematological screen MCV + MCH and HbA2		
	Sensitivity	Specificity
α-thalassemia	16.53%	97.24%
β-thalassemia	73.74%	96.10%

Information derived from reference [1]

Due to the complex nature of thalassemia genetic variants, conventional molecular analysis methods cannot provide a comprehensive solution, potentially leading to missed diagnoses, misdiagnoses, or delays in the diagnostic process.

Conventional PCR-based methods, such as PCR-RDB and Gap-PCR, typically focus on a limited range of major or common variants. Although next-generation sequencing (NGS) expands the coverage of variants, its ability to detect structural variants and rare deletions is still constrained by the technical limitations of short-read sequencing.

Comprehensive Analysis of Thalassemia Alleles (CATSA)

shapes the future of thalassemia screening

With single molecule real-time (SMRT) sequencing, a comprehensive solution for thalassemia genetic analysis has been developed — CATSA. This innovative approach simplifies the diagnostic process into a single solution, eliminating the need for multiple tests and significantly reducing diagnostic time.

Patent (China): ZL 202110329821.X

- Covers >2,200 variants to reliably identify rare carriers across high-diversity and isolated communities
- Determines *cis/trans* configuration to resolve phasing ambiguity, ensuring accurate zygosity reporting in multi-mutation cases

Technology: Single molecule real-time (SMRT) sequencing

Platform: PacBio Vega, Sequel II, and Sequel IIe system

Sample type: Dried blood spot (DBS), gDNA, blood, buccal swab, and amniotic fluid

Operation hour: 77 hours


(starting from sample preparation to report generation)

Expert consensus officially endorsed as a diagnostic and screening tool by the Medical Genetics Branch of the Chinese Medical Doctor Association (2025) [2]

	PCR-based assays				Sequencing	
	Gap-PCR	PCR-RDB	SMelting curve	+MLPA	NGS	CATSA
Deletions	✓	✓	✓	✓	✓	✓
Point Mutations	×	✓	✓	×	✓	✓
Homologous Recombinations	×	×	×	✓	×	✓

+ Applicable only for large deletion

§ Applicable only for common deletions & point mutations for α-thalassemia as well as common point mutations for β-thalassemia



Comprehensive Analysis of Thalassemia Alleles

Category	Basic panel	Expanded panel	Comprehensive panel	Prenatal comprehensive panel
α-thalassemia	30 Deletions	30 Deletions	30 Deletions	30 Deletions
	4 Homologous recombination	4 Homologous recombination	5 Homologous recombination	5 Homologous recombination
	80 <i>HBA1</i> SNVs/InDels; 139 <i>HBA2</i> SNVs/InDels	80 <i>HBA1</i> SNVs/InDels; 139 <i>HBA2</i> SNVs/InDels	80 <i>HBA1</i> SNVs/InDels; 139 <i>HBA2</i> SNVs/InDels	80 <i>HBA1</i> SNVs/InDels; 139 <i>HBA2</i> SNVs/InDels
β-thalassemia	28 Deletions	28 Deletions	28 Deletions	28 Deletions
	340 <i>HBB</i> SNVs/InDels	340 <i>HBB</i> SNVs/InDels	340 <i>HBB</i> SNVs/InDels	340 <i>HBB</i> SNVs/InDels
Hemoglobinopathy	--	324 <i>HBA1</i> SNVs/InDels; 386 <i>HBA2</i> SNVs/InDels	324 <i>HBA1</i> SNVs/InDels; 386 <i>HBA2</i> SNVs/InDels	324 <i>HBA1</i> SNVs/InDels; 386 <i>HBA2</i> SNVs/InDels
	--	694 <i>HBB</i> SNVs/InDels	694 <i>HBB</i> SNVs/InDels	694 <i>HBB</i> SNVs/InDels
Others	--	--	160 <i>HBD</i> SNVs/InDels	160 <i>HBD</i> SNVs/InDels
	--	--	102 SNVs/InDels in the <i>HBA1/2</i> and <i>HBB</i> genes that are low frequency, not clearly associated	<i>HBA1/2</i> and <i>HBB</i> variants that are low-frequency, not clearly associated

Performance backed by published data

The new gold standard

In a prospective multicenter study, CATSA demonstrated 100% accuracy in detecting pathogenic thalassemia variants among 1,159 individuals with abnormal hemoglobin parameters. [3] Additionally, in a separate study with 1,122 individuals, CATSA outperformed NGS by identifying and correcting discordant results with 100% accuracy, driving a 2.28% improvement in overall detection yield. [4]

Improved pregnancy outcomes

In 278 at-risk amniotic fluid samples, CATSA outperformed PCR-based assays by accurately identifying all genotypes, and correcting the predicted phenotype severity in 8 fetuses. Notably, 2 cases were reclassified from thalassemia trait to intermedia and one from intermedia to trait. [5]

Case study: A couple undergoes thalassemia carrier screening before giving birth

1. Post-test genetic counselling: low-risk and regular follow-up

	Mother	Father
Complete blood count	Anemia (Low MCH and MCV)	Normal
PCR test	β-thalassemia carrier (α / αα; β ^{CD71/72} / β ^N)	α-thalassemia carrier (-α ^{3.7} / αα; β ^N / β ^N)

2. However, their baby shows severe thalassemia symptoms at 8 months old.

- PCR results show baby genotype is concordant with parental genotypes
- Parents are also different types of thalassemia carrier

Question: Why is the baby exhibiting thalassemia symptoms when the genotype suggests he should not?

3. After having an advanced CATSA test, it is clear that the incorrect results from the father led to misdiagnosis.

	Father	Baby
PCR test	-α ^{3.7} / αα; β ^N / β ^N	-α ^{3.7} / αα; β ^{CD71/72} / β ^N
CATSA	-α ^{3.7} / αα ^{anti3.7} ; β ^N / β ^N	-α ^{3.7} / αα ^{anti3.7} ; β ^{CD71/72} / β ^N

This complex case of β-thalassemia heterogeneous (β^{CD71/72} / β^N) with α-triplicate (αα^{anti3.7}) enlarge the imbalance of α/β globin ratio, leading to β-thalassemia intermedia.

- The previous tests cannot detect α-triplicate (αα^{anti3.7})
- This couple has 25% chance of getting β-thalassemia intermedia offspring

Information in this document is subject to change without notice. Information provided is intended for reference only. All trademarks are the property of Berry Genomics, or their respective owners. References: [1] He J, Song W, Yang J, et al. Next-generation sequencing improves thalassemia carrier screening among premarital adults in a high prevalence population: the Dai nationality, China. *Genet Med*. 2017;19(9):1022-1031. [2] Consortium For The Application Of Single-Molecule Real-Time Sequencing For The Precision Medicine And Control Of Thalassemia, Group Of Clinical Genetics Medical Genetics Branch Of Chinese Medical Doctor Association, Wu L. *Zhonghua Yi Xue Yi Chuan Xue Za Zhi*. 2025;42(4):385-396. [3] Liang Q, Gu W, Chen P, et al. A More Universal Approach to Comprehensive Analysis of Thalassemia Alleles (CATSA). *J Mol Diagn*. 2021;23(9):1195-1204. [4] Huang R, Liu Y, Xu J, et al. Back-to-Back Comparison of Third-Generation Sequencing and Next-Generation Sequencing in Carrier Screening of Thalassemia. *Arch Pathol Lab Med*. 2024;148(7):797-804. [5] Liang Q, He J, Li Q, et al. Evaluating the Clinical Utility of a Long-Read Sequencing-Based Approach in Prenatal Diagnosis of Thalassemia. *Clin Chem*. 2023;69(3):239-250. INT_TGSthalav2022026