

Reaching further into dynamic mutations with the long-read era


# Tandem repeats

## Background

Tandem repeats are a common type of genetic variation in the human genome, ranging from 2 bp to over 100 bp. More than 60 tandem repeat loci have been linked to over 69 diseases [1], including Fragile X syndrome, Huntington's disease, and spinocerebellar ataxias.

Many of these conditions exhibit high phenotypic variability and overlapping clinical features, making diagnosis based on symptoms alone unreliable and sometimes leading to incorrect molecular testing approaches. A precise, multiplex molecular assay is needed, but routine tests face technological limitations, often leaving cases unresolved.

	Southern blot	RP-PCR	NGS
Extreme GC content	✓	✓	×
Interruptive motifs	×	×	×
Large repeat units (VNTR)	✓	×	×
Number of gene per assay	1	1-10	10-30
Repeat count accuracy	Low	Low (Limited to count <200)	Low (Limited to size <150 bp)
Throughput	Low	Low	High



## Challenge of routine assays:

### Interruptive motifs

Not designed to detect interruptive accurately

### Variability in repeat sizes (from STR to VNTR)

RP-PCR is limited to STRs (usually 2–6 bp)

NGS short-read length (~150 bp) is insufficient to span large expansion

### Flanking region mutations

Off-target may cause false results in PCR assays

### Extreme GC content of repeat units

Limiting the reliability of NGS because of low coverage

### Most located in non-coding regions

Blind spots for whole-exome sequencing (WES)

## Long-read sequencing opens new frontiers in the detection

To date, nearly half of the disorders associated with expansion loci have been discovered in the past decade, as previous methodological challenges hindered their detection.

Advances in long-read sequencing have accelerated the detection of expansions, eliminating the need for linkage data or prior hypotheses about the expanded motif.

## Timeline of repeat expansion discovery in human disorders




Figure adopted from reference [1] with modification

# dmTGS - massive tandem repeat panel with long-read sequencing

## dmTGS - massive tandem repeat panel with long-read sequencing

pioneers the future of dynamic mutation detection

dmTGS overcomes these challenges using HiFi reads with SMRT sequencing. This advanced approach surpasses PCR-based assays and short-read sequencing, providing a more comprehensive and reliable genetic profile for disease diagnosis.

**Technology:** Single molecule real-time (SMRT) sequencing  
**Platform:** PacBio Vega, Sequel II, and Sequel IIE system  
**Sample type:** Blood and gDNA  
**Turnaround time:** 29 working days  
 (starting from the date of sample arrival at the testing laboratory)

### HiFi reads

Accuracy > 99.9%  
 Long-read  
 Uniform coverage  
 No GC bias

### dmTGS

Repeat expansions  
 63 genes, 68 diseases

### dmTGS methodologies




Figure adopted from reference [2] with modification

### Comprehensive detection

Covers 68 diseases associated with 63 genes in a single assay, including the hard-to-map regions

### High resolution

Precisely distinguishes pathogenic from non-pathogenic repeat counts at a single-repeat level, reducing uncertainty in borderline cases

### Sensitive to low-level mosaicism

Detects mosaicism at low levels (e.g. 1% for *FMR1* CGG premutations and 5% for full mutations [2])

### Detailed characterization of polymorphic sequences

- Highly polymorphic sequences (e.g. VNTR) sometimes span dozens of kilobases
- Posing challenges for PCR-based assays (large repeat unit) and nanopore (high base error rates)
- dmTGS, with its single-nucleotide resolution, overcomes these challenges by enabling accurate mapping, repeat calling, and precise delineation of complex sequences and motifs

### Accurate detection of interruptive motifs


Detects interruptions within repetitive sequences that can impact genetic stability, transmission, and pathogenicity

For example:

- For Fragile X syndrome, women with AGG interruptions in CGG repeats are less likely to have children with full mutations [3]
- Similar mechanisms are observed in interruptive motifs in genes such as *ATXN1*, *ATXN2*, *DMPK*, and *HTT* [4, 5]


# dmTGS - massive tandem repeat panel with long-read sequencing

## Analysis of the 99-mer VNTR in *PLIN4* among 10 patients



PCR assays amplification can only provide an **approximate estimate of the repeat number** based on fragment size

In contrast, dmTGS accurately determined repeat counts in heterogeneous repeats. Sequence variability in the *PLIN4* gene is illustrated in composition plots.



## In one notable patient (P53) with phenotypes highly suspected with myotonic dystrophy:


1. **RP-PCR** assay for the *DMPK* gene (myotonic dystrophy type I) indicated CTG repeat count of 5/13 (**negative**)
2. Repeat analysis using "Expansion Hunter" on **WES** data reported a repeat count of 13/13 (**negative**)
3. **dmTGS** analysis detected a repeat count of 13/586 (**positive**) and identified a **large CCG interrupting motif** at the 3' end of the repeats, potentially **interfering with RP-PCR** extension

Information and figure derived from reference [2]

## Case study

### Clinical Information

- I-1 (Proband): Slurred speech, choking easily when drinking water, urinary incontinence, weakness in both lower limbs
- I-2 (Wife): Slurred speech, difficulty walking
- II-1 (Daughter): Asymptomatic
- II-2 (Elder son): Slurred speech, difficulty walking
- II-3 (Younger son): Slurred speech, difficulty walking



### Previous Genetic Test Results

1. **WES + CNV-seq** for I-1, I-2, and II-2: **Negative**
2. **RP-PCR** panel for spinocerebellar ataxia (10 genes) for I-1: Spinocerebellar ataxia 3 (SCA3) positive
3. SCA3 test for other members:
  - o II-1: Negative
  - o II-2: Positive
  - o II-3: **Negative**

The underlying cause of the phenotype in I-2 and II-3 remains unknown.

### To further investigate the cause, dmTGS is performed:

- I-2 and II-3: Neuronal intranuclear inclusion disease (NIID) (*NOTCH2NLC* 22/68) positive
- II-2: SCA3 (*ATXN3* 16/73) positive and neuronal intranuclear inclusion disease (NIID) (*NOTCH2NLC* 22/69) positive
- Abnormal expansion in the *NOTCH2NLC* gene is associated with NIID, which matches the clinical phenotypes of the I-2 and II-3

In the Chinese population, SCA3 is the most common ataxia. While SCA3 or SCA10 panel testing is preferred for typical cases, overlapping phenotypes in tandem repeat disorders make symptom-based diagnosis unreliable.

**dmTGS provides a comprehensive single-test solution, minimizing misdiagnosis risk.**

References: [1] Depienne C, Mandel JL. 30 years of repeat expansion disorders: What have we learned and what are the remaining challenges?. *Am J Hum Genet.* 2021;108(5):764-785. [2] Yang K, Liu Y, Zhang J, et al. dmTGS: Precise Targeted Enrichment Long-Read Sequencing Panel for Tandem Repeat Detection. *Clin Chem.* 2025;71(2):319-331. [3] Yrigollen CM, Durbin-Johnson B, Gane L, Nelson DL, Hagerman R, Hagerman PJ, Tassone F. AGG interruptions within the maternal *FMR1* gene reduce the risk of offspring with fragile X syndrome. *Genet Med.* 2012 Aug;14(8):729-36. [4] Menon RP, Nethisinghe S, Faggiano S, et al. The role of interruptions in polyQ in the pathology of SCA1. *PLoS Genet.* 2013;9(7):e1003648. [5] Nethisinghe S, Pigazzini ML, Pemble S, et al. PolyQ Tract Toxicity in SCA1 is Length Dependent in the Absence of CAG Repeat Interruption. *Front Cell Neurosci.* 2018;12:200.

# dmTGS - massive tandem repeat panel with long-read sequencing

## List of disorders covered by dmTGS

### Ataxia (20)

Diseases	Gene	Repeat Unit
Spinocerebellar ataxia 1 (SCA1)	ATXN1	CAG, CAT
Spinocerebellar ataxia 2 (SCA2)	ATXN2	CAG, CAA
Spinocerebellar ataxia 3 (SCA3)	ATXN3	CAG, CAA
Spinocerebellar ataxia 4 (SCA4)	ZFH3	GGC
Spinocerebellar ataxia 6 (SCA6)	CACNA1A	CAG
Spinocerebellar ataxia 7 (SCA7)	ATXN7	CAG
Spinocerebellar ataxia 8 (SCA8)	ATXN8OS_ATXN8	CTG-CAG
Spinocerebellar ataxia 10 (SCA10)	ATXN10	ATTCT
Spinocerebellar ataxia 12 (SCA12)	PPP2R2B	CAG
Spinocerebellar ataxia 17 (SCA17)	TBP	CAG, CAA
Spinocerebellar ataxia 31 (SCA31)	BEAN1, TK2	TGGAA
Spinocerebellar ataxia 27B, late-onset (SCA27B)	FGF14	GAA
Spinocerebellar ataxia 36 (SCA36)	NOP56	GGCCTG
Spinocerebellar ataxia 37 (SCA37)	DAB1	ATTTT
Spinocerebellar ataxia 51 (SCA51)	THAP11	CAG
Dentatorubral-pallidolusian atrophy (DRPLA)	ATN1	CAG
Friedreich ataxia (FRDA)	FXN	GAA
Global developmental delay, progressive ataxia, and elevated glutamine (GDPAG)	GLS	GCA
Cerebellar ataxia, neuropathy, and vestibular areflexia syndrome (CANVAS)	RFC1	AAGGG
Fragile X tremor/ataxia syndrome (FXTAS)	FMR1	CGG, AGG

### Epilepsy (8)

Diseases	Gene	Repeat Unit
Epilepsy, progressive myoclonic 1A (EPM1A)	CSTB	CCCCGCCCGCGC
Epilepsy, familial adult myoclonic, 1 (FAME1)	SAMD12	TTTCA/TTTTA
Epilepsy, familial adult myoclonic, 2 (FAME2)	STARD7	AAATG
Epilepsy, familial adult myoclonic, 3 (FAME3)	MARCHF6	TTTCA/TTTTA
Epilepsy, familial adult myoclonic, 4 (FAME4)	YEATS2	TTTCA
Epilepsy, familial adult myoclonic, 6 (FAME6)	TNRC6A	TTTCA
Epilepsy, familial adult myoclonic, 7 (FAME7)	RAPGEF2	TTTCA
Tremor, hereditary essential, 6 (ETM6)	NOTCH2NLC	GGC

### Developmental disorders (8)

Diseases	Gene	Repeat Unit
Synpolydactyly 1 (SPD1)	HOXD13	GCN
Richieri-Costa-Pereira syndrome (RCPS)	EIF4A3	18-20 bp
Hand-foot-genital syndrome (HFGS)	HOXA13	GCN
Blepharophimosis, ptosis, and epicanthus inversus (BPES)	FOXL2	GCN
VACTERL association, X-linked	ZIC3	GCN
Oculo-auriculo-vertebral spectrum (OVAS)		
Tetralogy of Fallot (ToF)	TBX1	GCN
Baratela-Scott syndrome (BSS)	XYLT1	GGC

For details, please contact [marketing@xcelom.com](mailto:marketing@xcelom.com)

Information in this document is subject to change without notice. Information provided is intended for reference only. All trademarks are the property of Berry Genomics, or their respective owners.

INT\_TGSdmTGS2v2.02012026

### Other neurological disorders (18)

Diseases	Gene	Repeat Unit
Developmental and epileptic encephalopathy 1 (DEE1)	ARX(ARX_1)	GCN
Intellectual developmental disorder, X-linked 29 (XLID29)		
Fragile X syndrome (FXS)/ Fragile X-associated primary ovarian insufficiency (FXPOI)/ Fragile X-associated neuropsychiatric disorders (FXAND)	FMR1	CGG, AGG
Intellectual developmental disorder, X-linked 109 (XLID109)	AFF2	CCG
KINSHIP syndrome (KINS)	AFF3	CGG
Autism spectrum disorder (FRA7A)	ZNF713	CGG
Intellectual developmental disorder, autosomal dominant, FRA12A type (MRFRA12A)	DIP2B	CGG
Jacobsen syndrome (JBS)	CBL	CCG
Huntington disease (HD)	HTT	CAG
Huntington disease-like 2 (HDL2)	JPH3	CTG
Neuropathy, hereditary sensory and autonomic, type VIII (HSAN8)	PRDM12	GCC
Creutzfeldt-Jakob disease (CJD)	PRNP	24 bp
Intellectual developmental disorder, X-linked, with isolated growth hormone deficiency	SOX3	GCN
Panhypopituitarism, X-linked		
Holoprosencephaly 5 (HP5)	ZIC2	GCN
Developmental and epileptic encephalopathy 1 (DEE1)	ARX(ARX_2)	GCN
Partington syndrome (PRTS)		
Intellectual developmental disorder, X-linked 29 (XLID29)		

### Neuromuscular disorders (12)

Diseases	Gene	Repeat Unit
Spinal and bulbar muscular atrophy, X-linked 1 (SBMA)	AR	CAG
Myotonic dystrophy 1 (DM1)	DMPK	CTG
Myotonic dystrophy 2 (DM2)	CNBP	CCTG, GCTC, TCTG
Oculopharyngodistal myopathy 3 (OPDM3)	NOTCH2NLC	GGC
Oculopharyngodistal myopathy 1 (OPDM1)	LRP12	CGG
Oculopharyngodistal myopathy 2 (OPDM2)	GIPC1	CGG
Oculopharyngodistal myopathy 4 (OPDM4)	RILPL1	CGG
Oculopharyngeal myopathy with leukoencephalopathy 1 (OPML1)	NUTM2BAS1 LOC642361	CCG-CGG
Oculopharyngeal muscular dystrophy (OPMD)	PABPN1	GCN
Neuronal intranuclear inclusion disease (NIID)	NOTCH2NLC	GGC
Frontotemporal dementia and/or amyotrophic lateral sclerosis 1 (FTDALS1)	C9orf72	GGGGCC
Myopathy with rimmed ubiquitin-positive autophagic vacuolation, autosomal dominant (MRUPAV)	PLIN4	99 bp

### Ophthalmological disorder (1)

Diseases	Gene	Repeat Unit
Corneal dystrophy, Fuchs endothelial, 3 (FECD3)	TCF4	CTG

### Respiratory disorder (1)

Diseases	Gene	Repeat Unit
Central hypoventilation syndrome, congenital, 1, with or without Hirschsprung disease (CCHS1)	PHOX2B	GCN